


Genetic Technology Instrumentation

GTI-Protein Electrophoresis Tank

Electrophoresis is used for the separation of charged molecules like proteins by an electric current applied to a matrix. Electrophoresis of haemoglobin and serum proteins can provide important diagnostic information. GTI Protein electrophoresis tank is specially designed to accommodate cellulose acetate membrane strips.

- Ideal for Cellulose Acetate Membrane (CAM) electrophoresis of haemoglobin and serum/plasma proteins
- Accommodates 6 X 15 cm CAM strips



**Genetic Technology Instrumentation
(GTI)**

Phone: +92-51-5167312

E-mail: info@grcpk.com

www.grcpk.com